




PRACTICAL 5 BLAST


唐凯临
2021.4




Score Matrix

```
EBLOSUM30
EBLOSUM35
EBLOSUM40
EBLOSUM45
EBLOSUM50
EBLOSUM55
EBLOSUM60
EBLOSUM62
EBLOSUM62-12
EBLOSUM65
EBLOSUM70
EBLOSUM75
EBLOSUM80
EBLOSUM85
EBLOSUM90
EBLOSUMN
EPAM10
EPAM100
EPAM110
EPAM120
EPAM130
EPAM140
EPAM150
EPAM160
EPAM170
EPAM180
EPAM190
EPAM20
EPAM200
EPAM210
EPAM220
EPAM230
EPAM240
EPAM250
EPAM260
EPAM270
EPAM280
EPAM290
EPAM30
EPAM300
EPAM310
EPAM320
EPAM330
EPAM340
EPAM350
EPAM360
EPAM370
EPAM380
EPAM390
EPAM40
EPAM400
EPAM410
EPAM420
EPAM430
EPAM440
EPAM450
EPAM460
EPAM470
EPAM480
EPAM490
EPAM50
EPAM500
EPAM60
EPAM70
EPAM80
EPAM90
```


Distribution of the top 310 Blast Hits on 100 subject sequences




Distribution of the top 285 Blast Hits on 100 subject sequences



Distribution of the top 341 Blast Hits on 100 subject sequences



Distribution of the top 31 Blast Hits on 30 subject sequence



REVIEW

Query Sequence	Search Database	BLAST Program	Sequence Comparison	BLAST output
DNA	nucleotide	blastn	compare query nucleotide against nucleotide db	Nucleotide
DNA	protein	blastx	translate query seq in all reading frames into amino acids, then compare with protein db	Amino acid
DNA	nucleotide	tblastx	translate both query & db seq in all reading frames, then compare between protein seqs	Amino acid
Protein	protein	blastp	compare query protein against protein db	Amino acid
Protein	nucleotide	tblastn	translate db nucleotide seq in all reading frames, then compare between protein seqs	Amino acid



PREVIEW

- Which is better if you want to find representative sequences from many species?
- Which is better if you want to find all related sequences from a single species?



REVIEW



- **nr/nt** database contains **ALL known sequences reported at NCBI**
- NCBI created two databases called **RefSeq_Protein** and **RefSeq_Genomic**, designed to **reduce duplication in nr/nt** by selecting unique representative sequences for each locus
- **Swissprot or Uniprot** is a database of **highly curated protein sequences** , representing an effort to annotate/enrich all the protein sequence records in **nr**



CUSTOMIZING BLAST TO YOUR NEEDS


- Build your own searchable database for a customized dataset and perform BLAST search against it
 - a set of your own sequences that are not available in the public databases.
 - eg, novel sequences from sequencing projects by your lab or collaborators.
 - a set of sequences that are available in public databases, but have been processed/organized/grouped according to your liking.
 - eg, protein sequences grouped according to their function.

DOS——DISK OPERATION SYSTEM



常用命令

- `dir [C:][path][filename][.ext][/o][/s][/p][/w][/a]`
- `md [C:]path`
- `cd [C:][path]`, `cd ..` 返回上层目录
- `rd [d:]path` 不能删除非空目录；不能删除当前目录。



```
Microsoft(R) Windows Millennium
(C)Copyright Microsoft Corp 1981-1999.

C:\WINDOWS\Desktop> md purple —— 创建一个新目录
C:\WINDOWS\Desktop> cd purple —— 进入建立的新目录
C:\WINDOWS\Desktop\purple> cd.. —— 退出Purple目录
C:\WINDOWS\Desktop> rd purple —— 删除Purple目录
C:\WINDOWS\Desktop> cd purple —— 再次进入，目录已被删除！
Invalid directory
C:\WINDOWS\Desktop>
```

常用命令

- `copy [C:][path][filename.ext] [C:][path]filename.ext`
- `del [C:][path]filename.ext`
- `ren [C:][path]filename1[.ext] filename2[.ext]`
- `cls` 清除
- 系统命令 + `/?` 帮助
- 通配符: * 和 ?
 - * 表示一个字符串
 - ? 只代表一个字符



如何本地运行BLAST

- ftp://ftp.ncbi.nlm.nih.gov/blast/



DOWNLOAD



FTP 目录 /blast/executables/blast+/LATEST 位于 f

若要在文件资源管理器中查看此 FTP 站点，请单击“视图”，然后单击“在文件资源

[转到高层目录](#)

12/04/2019 02:52上午	85	ChangeLog
12/04/2019 02:50上午	20,367,036	ncbi-blast-2.10.0+-4.src.rpm
12/04/2019 02:53上午	63	ncbi-blast-2.10.0+-4.src.rpm.md5
12/04/2019 02:50上午	183,553,344	ncbi-blast-2.10.0+-4.x86_64.rpm
12/04/2019 02:53上午	66	ncbi-blast-2.10.0+-4.x86_64.rpm.md5
12/04/2019 02:53上午	25,547,460	ncbi-blast-2.10.0+-src.tar.gz
12/04/2019 02:53上午	64	ncbi-blast-2.10.0+-src.tar.gz.md5
12/04/2019 02:53上午	29,920,102	ncbi-blast-2.10.0+-src.zip
12/04/2019 02:53上午	61	ncbi-blast-2.10.0+-src.zip.md5
12/04/2019 02:49上午	90,788,089	ncbi-blast-2.10.0+-win64.exe
12/04/2019 02:53上午	63	ncbi-blast-2.10.0+-win64.exe.md5
12/04/2019 02:52上午	233,258,021	ncbi-blast-2.10.0+-x64-linux.tar.gz
12/04/2019 02:53上午	70	ncbi-blast-2.10.0+-x64-linux.tar.gz.md5
12/04/2019 02:53上午	147,458,501	ncbi-blast-2.10.0+-x64-macosx.tar.gz
12/04/2019 02:53上午	71	ncbi-blast-2.10.0+-x64-macosx.tar.gz.md5
12/04/2019 02:50上午	90,505,163	ncbi-blast-2.10.0+-x64-win64.tar.gz
12/04/2019 02:53上午	70	ncbi-blast-2.10.0+-x64-win64.tar.gz.md5
12/04/2019 02:52上午	149,443,790	ncbi-blast-2.10.0+.dmg
12/04/2019 02:53上午	57	ncbi-blast-2.10.0+.dmg.md5

DOWNLOAD

○ <ftp://ftp.ncbi.nlm.nih.gov/blast/db/>

04/23/2020	02:51下午	2, 235, 750, 412	refseq protein.11.tar.gz
04/23/2020	02:51下午	59	refseq protein.11.tar.gz.md5
04/23/2020	02:51下午	2, 235, 686, 988	refseq protein.12.tar.gz
04/23/2020	02:51下午	59	refseq protein.12.tar.gz.md5
04/23/2020	02:52下午	2, 235, 708, 007	refseq protein.13.tar.gz
04/23/2020	02:52下午	59	refseq protein.13.tar.gz.md5
04/23/2020	02:52下午	2, 235, 784, 183	refseq protein.14.tar.gz
04/23/2020	02:52下午	59	refseq protein.14.tar.gz.md5
04/23/2020	02:52下午	2, 235, 733, 660	refseq protein.15.tar.gz
04/23/2020	02:52下午	59	refseq protein.15.tar.gz.md5
04/23/2020	02:53下午	2, 235, 730, 728	refseq protein.16.tar.gz
04/23/2020	02:53下午	59	refseq protein.16.tar.gz.md5
04/23/2020	02:53下午	2, 235, 679, 301	refseq protein.17.tar.gz
04/23/2020	02:53下午	59	refseq protein.17.tar.gz.md5
04/23/2020	02:53下午	2, 235, 391, 018	refseq protein.18.tar.gz
04/23/2020	02:53下午	59	refseq protein.18.tar.gz.md5
04/23/2020	02:54下午	2, 235, 879, 501	refseq protein.19.tar.gz
04/23/2020	02:54下午	59	refseq protein.19.tar.gz.md5
04/23/2020	02:54下午	2, 198, 649, 196	refseq protein.20.tar.gz
04/23/2020	02:54下午	59	refseq protein.20.tar.gz.md5
04/24/2020	05:03下午	3, 214, 671, 722	refseq rna.00.tar.gz
04/24/2020	05:03下午	55	refseq rna.00.tar.gz.md5
04/24/2020	05:03下午	2, 400, 117, 599	refseq rna.01.tar.gz
04/24/2020	05:03下午	55	refseq rna.01.tar.gz.md5
04/24/2020	05:04下午	2, 336, 666, 179	refseq rna.02.tar.gz
04/24/2020	05:04下午	55	refseq rna.02.tar.gz.md5
04/24/2020	05:04下午	2, 269, 547, 740	refseq rna.03.tar.gz
04/24/2020	05:04下午	55	refseq rna.03.tar.gz.md5
04/24/2020	05:05下午	2, 227, 633, 648	refseq rna.04.tar.gz
04/24/2020	05:05下午	55	refseq rna.04.tar.gz.md5


RUN LOCAL BLAST CLIENT

- Download the correct version of BLAST for your computer
- Install BLAST;
- copy query & database fasta files into your folder
- Customize BLAST database & run BLAST search through command window

Name	Date modified	Type	Size	Tags
blast_formatter.exe	14/8/2010 2:25 AM	Application	6,264 KB	
blastdb.fasta	11/9/2010 3:35 AM	FASTA File	1 KB	
blastdb_aliastool.exe	14/8/2010 2:25 AM	Application	1,804 KB	
blastdbcheck.exe	14/8/2010 2:25 AM	Application	2,868 KB	
blastdbcmd.exe	14/8/2010 2:25 AM	Application	4,104 KB	
blastn.exe	14/8/2010 2:25 AM	Application	6,388 KB	
blastp.exe	14/8/2010 2:25 AM	Application	6,384 KB	
blastx.exe	14/8/2010 2:25 AM	Application	6,372 KB	

改变路径

- 改变目录: cd
- 改变盘符

 命令提示符

```
C:\Users\tangk>cd D:\Program Files\NCBI\blast-2.2.31+\bin  
C:\Users\tangk>d:  
D:\Program Files\NCBI\blast-2.2.31+\bin>
```


COMMAND

○makeblastdb -help

```
D:\test\blast-2.2.30+\bin>makeblastdb -help
USAGE
  makeblastdb.exe [-h] [-help] [-in input_file] [-input_type type]
                  -dbtype molecule_type [-title database_title] [-parse_seqids]
                  [-hash_index] [-mask_data mask_data_files] [-mask_id mask_algo_ids]
                  [-mask_desc mask_algo_descriptions] [-gi_mask]
                  [-gi_mask_name gi_based_mask_names] [-out database_name]
                  [-max_file_sz number_of_bytes] [-logfile File_Name] [-taxid TaxID]
                  [-taxid_map TaxIDMapFile] [-version]

DESCRIPTION
  Application to create BLAST databases, version 2.2.30+

REQUIRED ARGUMENTS
  -dbtype <String, 'nucl', 'prot'>
    Molecule type of target db

OPTIONAL ARGUMENTS
  -h
    Print USAGE and DESCRIPTION; ignore all other parameters
  -help
    Print USAGE, DESCRIPTION and ARGUMENTS; ignore all other parameters
```

COMMAND

blastp -help

```
D:\test\blast-2.2.30+\bin>blastp -help
USAGE
blastp [-h] [-help] [-import_search_strategy filename]
        [-export_search_strategy filename] [-task task_name] [-db database_name]
        [-dbsize num_letters] [-gilist filename] [-seqidlist filename]
        [-negative_gilist filename] [-entrez_query entrez_query]
        [-db_soft_mask filtering_algorithm] [-db_hard_mask filtering_algorithm]
        [-subject subject_input_file] [-subject_loc range] [-query input_file]
        [-out output_file] [-evalue evalue] [-word_size int_value]
        [-gapopen open_penalty] [-gapextend extend_penalty]
        [-qcov_hsp_perc float_value] [-xdrop_ungap float_value]
        [-xdrop_gap float_value] [-xdrop_gap_final float_value]
        [-searchsp int_value] [-max_hsps int_value] [-sum_stats bool_value]
        [-seg SEG_options] [-soft_masking soft_masking] [-matrix matrix_name]
        [-threshold float_value] [-culling_limit int_value]
        [-best_hit_overhang float_value] [-best_hit_score_edge float_value]
        [-window_size int_value] [-lcase_masking] [-query_loc range]
        [-parse_deflines] [-outfmt format] [-show_gis]
        [-num_descriptions int_value] [-num_alignments int_value]
        [-line_length line_length] [-html] [-max_target_seqs num_sequences]
        [-num_threads int_value] [-ungapped] [-remote] [-comp_based_stats compo]
```



建库

○makeblastdb -in test\ecoli.aa -dbtype prot -out test\ecolidb.fasta

```
D:\Program Files\NCBI\blast-2.2.31+\bin>makeblastdb -in test\ecoli.aa -dbtype prot -out test\ecolidb.fasta

Building a new DB, current time: 03/12/2020 20:46:20
New DB name:   D:\Program Files\NCBI\blast-2.2.31+\bin\test\ecolidb.fasta
New DB title:  test\ecoli.aa
Sequence type: Protein
Keep Linkouts: T
Keep MBits: T
Maximum file size: 1000000000B
Adding sequences from FASTA; added 4289 sequences in 0.140971 seconds.
```



比对

o `blastp -db test\ecolidb.fasta -query test\myecoliquery.txt -num_alignments 1 -evaluate 1e-5 -out test\ecoliout`

```
D:\Program Files\NCBI\blast-2.2.31+\bin>blastp -db test\ecolidb.fasta -query test\myecoliquery.txt -evaluate 1e-5 -out test\ecoliout  
D:\Program Files\NCBI\blast-2.2.31+\bin>
```

DATA (D:) > Program Files > NCBI > blast-2.2.31+ > bin > test

名称 ^	
★	ecoli.aa
★	ecolidb.fasta.phr
★	ecolidb.fasta.pin
★	ecolidb.fasta.psq
★	ecoliout
	myecoliquery

结果

BLASTP 2.2.31+

Reference: Stephen F. Altschul, Thomas L. Madden, Alejandro A. Schaffer, Jinghui Zhang, Zheng Zhang, Webb Miller, and David J. Lipman (1997), "Gapped BLAST and PSI-BLAST: a new generation of protein database search programs", Nucleic Acids Res. 25:3389-3402.

Reference for composition-based statistics: Alejandro A. Schaffer, L. Aravind, Thomas L. Madden, Sergei Shavirin, John L. Spouge, Yuri I. Wolf, Eugene V. Koonin, and Stephen F. Altschul (2001), "Improving the accuracy of PSI-BLAST protein database searches with composition-based statistics and other refinements", Nucleic Acids Res. 29:2994-3005.

Database: test\ecoli.aa
4,289 sequences; 1,358,990 total letters

Query= My_peptide

Length=640

	Score	E
Sequences producing significant alignments:	(Bits)	Value
gi 1786539 gb AAC73447.1 (AE000141) beta-D-galactosidase [Esch...	1329	0.0
gi 1789457 gb AAC76111.1 (AE000389) evolved beta-D-galactosida...	406	5e-129
gi 1787903 gb AAC74689.1 (AE000257) beta-D-glucuronidase [Esch...	101	2e-023



TIPS

○ Too much results

- Refseq
- Weight matrix
- $E \downarrow$
-

○ Less results

- Nr/nt
- $E \uparrow$
- Weight matrix
- Word \downarrow
-



基本的BLAST与特别的BLAST

nucleotide blast	Search a nucleotide database using a nucleotide query <i>Algorithms: blastn, megablast, discontinuous megablast</i>
protein blast	Search protein database using a protein query <i>Algorithms: blastp, psi-blast, phi-blast, delta-blast</i>
blastx	Search protein database using a translated nucleotide query
tblastn	Search translated nucleotide database using a protein query
tblastx	Search translated nucleotide database using a translated nucleotide query

Specialized BLAST

Choose a type of specialized search (or database name in parentheses.)

- ❑ Make specific primers with [Primer-BLAST](#)
- ❑ Cluster multiple sequences together with their database neighbors using [MOLE-BLAST](#)
- ❑ Find [conserved domains](#) in your sequence (cds)
- ❑ Find sequences with similar [conserved domain architecture](#) (cdart)
- ❑ Search sequences that have [gene expression profiles](#) (GEO)
- ❑ Search [immunoglobulins and T cell receptor sequences](#) (IgBLAST)
- ❑ Screen sequence for [vector contamination](#) (vecscreen)
- ❑ [Align](#) two (or more) sequences using BLAST (bl2seq)
- ❑ Search [protein](#) or [nucleotide](#) targets in PubChem BioAssay
- ❑ Search [SRA by experiment](#)
- ❑ Constraint Based Protein [Multiple Alignment Tool](#)
- ❑ Needleman-Wunsch [Global Sequence Alignment Tool](#)
- ❑ Search [RefSeqGene](#)
- ❑ Search [trace archives](#)
- ❑ Search bacterial and fungal rRNA sequences with [Targeted Loci BLAST](#)



SUMMARY OF BLAST

➤ Introduction on BLAST

- What is BLAST ?
- BLAST flavors ?
- BLAST databases ?
- BLAST Access: Web Blast ? Local BLAST Client?

➤ Application of BLAST

- How BLAST can be applied for life science research ?
 - eg, sequence identity verification; conserved domains; similar genes/proteins; distant relatives; homology; etc.
- Advantages/Disadvantages of BLAST