


Protein purification

王春光

chunguangwang@tongji.edu.cn



2020-09-29

myoglobin vs. hemoglobin



myoglobin: O₂ storage

hemoglobin:
O₂ transportation





Mb: myoglobin (肌红蛋白);
Hb: hemoglobin (血红蛋白).


Something else to know about hemoglobin

Hemoglobin binds CO with higher affinity (亲和力);

Sickle-cell anemia (镰刀形贫血症) caused by hemoglobin mutation $E\beta$




ATP synthase is a highly complex molecular machine.





(b)

Copyright 1999 John Wiley and Sons, Inc. All rights reserved.

ATP synthase is a highly complex molecular machine.





PDB



Basic principle of protein purification




The unique property matters!

Learning objectives:

- Explain the ion exchange chromatography（离子交换层析）；
- Explain the gel filtration chromatography（凝胶过滤层析）；
- Choose an appropriate purification technique.


Ion exchange chromatography 离子交换层析

The net charge of a protein as a function of pH




(from GE Healthcare)

Ion exchange chromatography 离子交换层析



Ion exchange chromatography 离子交换层析





Type of Ion Exchange	Common Abbreviation	Functional group	
Strong Anion	Q	Quarternary Ammonium	$-\text{CH}_2\text{N}^+(\text{CH}_3)_3$
Weak Anion	DEAE	Diethylaminoethyl	$-\text{C}_2\text{H}_4\text{N}^+\text{H}(\text{C}_2\text{H}_5)_2$
Strong Cation	S	Sulfonic Acid	$-\text{CH}_2\text{SO}_3^-$
Weak Cation	CM	Carboxymethyl	$-\text{CH}_2\text{COO}^-$

Elution (洗脱) of ion exchange chromatography

By increasing salt concentration or changing pH condition;

Stepwise elution (阶段洗脱) vs. gradient elution (梯度洗脱);



Example: Mono S column from GE Healthcare

Column Mono S 5/50 GL

Sample: 1. Ribonuclease A (1.5 mg/ml) SIGMA R-5000
2. Cytochrome C (0.4 mg/ml) SIGMA C-7752
3. Lysozyme (0.4 mg/ml) SIGMA L-6876


Sample volume: 100 μ l

Gradient: 0–100% elution buffer in 20 CV

Start buffer: 20 mM sodium phosphate, pH 6.8

Elution buffer: 20 mM sodium phosphate + 0.4 M NaCl, pH 6.8

Flow rate: 1.0 ml/min (room temperature)






Ribonuclease A
pI=9.45

Cytochrome C
pI=10.0-10.5


Lysozyme
pI=11.35

Hydrophobic interaction chromatography (HIC) (疏水相互作用层析)

Hydrophobicity: octyl > phenyl > butyl



(http://www.chisso.co.jp/fine/en/cellufine/images/hyd_partial.gif)




Lab Pack: Phenyl Sepharose High Performance
Column: 5 × 50 mm High Performance column
Bed height: 10 mm
Sample: 60 µg active and 60 µg methylamine treated, inactive α₂-macroglobulin
Start buffer (A): 0.02 M sodium phosphate with 1 M Na₂SO₄, pH 7.2
Elution buffer (B): 0.025 M sodium phosphate, pH 7.2
Flow rate: 0.5 ml/min

(from GE Healthcare)


Protein size

amino acid



Mr = 133

residue





Mr = 115

Average molecular weight of residues = 110.

Example: a protein of 22 KD → about 200 residues;

→ cDNA of about 600 bp.




Gel filtration chromatography 凝胶过滤层析





Scanning electron micrograph of an agarose gel. Magnification x 50,000. (Ref. Anders S. Medin, PhD Thesis, Uppsala University 1995)

东 华 门 地 区

- 1北梅竹胡同
- 2景丰胡同
- 3骑河楼北巷
- 4福禄巷
- 5箭杆胡同
- 6北池子三条
- 7北池子二条
- 8北池子头条
- 9文书馆巷
- 10万庆巷
- 11智德前巷
- 12普渡寺东巷
- 13普渡寺后巷
- 14普渡寺前巷
- 15灯笼库胡同
- 16磁器库北巷
- 17磁器库南巷
- 18磁器库胡同
- 19北库司胡同
- 20飞龙桥胡同
- 21东厂北巷
- 22安居里
- 23大草厂胡同
- 24南口袋胡同
- 25桂花胡同
- 26玉石胡同
- 27沙滩南巷
- 28敦厚里
- 29东胜胡同
- 30东单小二条
- 31官房大院
- 32小草厂胡同
- 33苏州胡同
- 34小苏州胡同
- 35合基厂三条
- 36北口袋胡同
- 37北官场胡同



Gel filtration chromatography 凝胶过滤层析



Scanning electron micrograph of an agarose gel. Magnification x 50,000. (Ref. Anders S. Medin, PhD Thesis, Uppsala University 1995)

Also known as size exclusion chromatography;

Exclusion limit (排阻限度) ;

Effective separation range (globular proteins);
(有效分离范围)

Column volume and void volume (外水体积) ;


Limited sample volume.

Example: Superdex 75 from GE Healthcare


Bed dimensions	10 × 300–310 mm	
Bed volume	Approximately 24 ml	
	Superdex 75	Superdex 200
Exclusion limit, M_r , globular proteins	Approx. 1×10^5	Approx. 1.3×10^6
Optimum separation range globular proteins, M_r	3 000–70 000	10 000–600 000
dextran	500–30 000	1000–100 000
Flow rate (water at room temperature) recommended	0.5–1.0 ml/min	0.25–0.75 ml/min

Column Superdex 75 10/300 GL

Sample:	1. BSA (M_r 67 000) 8 mg/ml
	2. Ovalbumin (M_r 43 000) 2.5 mg/ml
	3. Ribonuclease A (M_r 13 700) 5 mg/ml
	4. Aprotinin (M_r 6 512) 2 mg/ml
	5. Vitamin B12 (M_r 1355) 0.1 mg/ml
Sample volume:	500 μ l
Eluent:	0.05 M phosphate buffer, 0.15 M NaCl, pH 7.0
Flow rate:	0.4 ml/min, room temperature
Detection:	280 nm




Affinity chromatography 亲和层析



Flash

Immobilized metal chelating chromatography (IMAC)

(金属螯合层析)



More options of some other metal ions (Co^{2+} , Cu^{2+});

Popularly used for $(\text{His})_6$ -tagged proteins.

protein purification techniques

➤ Based on physical and chemical characteristics:

size gel filtration chromatography

charge ion exchange chromatography


hydrophobicity hydrophobic interaction chromatography

➤ Based on biological characteristics:


affinity chromatography




Gel filtration



Hydrophobic interaction



Ion exchange



Affinity

Change
condition
to elute

No

Yes

Yes

Yes

Sample
volume
limited

Yes

No


No

No

AKTA purifier 10 from GE Healthcare



AKTA purifier 10 from GE Healthcare



Next time:

Protein characterization.